

D1008 series Palm Micro Centrifuge

Features

- Ideal for microfiltration and quick spin-downs
- Easy-to-use with dual start/stop function
- 10000rpm fixed speed (D1008U)
- 7000rpm fixed speed (D1008)
- 5000rpm fixed. speed (D1008E)
- Compatible with 2mL x 8 and PCR8 x 2 rotors
- Quiet running, noise level ≤ 45 dB
- Rotor change without tools
- Accelerates and brakes in seconds
- Four colors available, green, blue, yellow and pink, add color to your routine laboratory life



Quiet operation



Vibration free



CE FC

Specifications

	D1008U	D1008	D1008E
Max. Speed	10000rpm	7000rpm	5000rpm
Max. RCF	5460×g	2680×g	1360×g
Rotor Capacity	0.2/0.5/1.5/2.0mL×8; 0.2mL×16 PCR strips or 0.2mL×2 PCR 8 strips	0.2/0.5/1.5/2.0mL×8; 0.2mL×16 PCR strips or 0.2mL×2 PCR 8 strips	0.2/0.5/1.5/2.0mL×8; 0.2mL×16 PCR strips or 0.2mL×2 PCR 8 strips
Run Time	Continuous operation	Continuous operation	Continuous operation
Driving Motor	DC motor	DC motor	DC motor
Power	AC110-240V/50Hz/60Hz	AC110-240V/50Hz/60Hz	AC110-240V/50Hz/60Hz
Noise level	≤ 45 dB	≤ 45 dB	≤ 45 dB
Dimensions [D×W×H]	160 x 170 x 122mm	160 x 170 x 122mm	160 x 170 x 122mm
Weight	0.5kg	0.5kg	0.5kg
Certification	CE cTUVus FCC MCA	CE cTUVus FCC MCA	CE cTUVus FCC MCA