

dTrite Electronic Titrator

The dTrite Digital Burette delivers accurate, precise and convenient bottle-top titration, as well as optimum operator safety.

It combines the magnetic stirrer with the titration process, making it easier for the operator to achieve the satisfactory results in the field of analytical chemistry, food industry, water analysis,, etc.



Features

- High quality motor providing excellent precision and accuracy
- Electronic control decreases the repetitive strain injury
- The remote control panel makes the operating process without the manual disturbance
- Build-in magnetic stirrer support

Specifications

Volume range	0.01mL-99.99mL, Max piston lift is 10mL, resolution is 10 μ L
Volume accuracy	R=0.2% CV=0.07%
Velocity	16 stages
Operating temperature range	10°C~30°C
Quality control	DIN EN ISO 8655
Control type	External control, including the stirrer and the burette program