

# Pipetting Filler

## Levo Plus Pipette Filler

### Features

- Appropriate for most of the plastic and glass pipettes from 0.1 -100ml
- Eight speeds of selection for aspiration and dispensing
- Large LCD display showing low battery warning and speed settings
- Enabling single-handed operation with minimum effort
- Light and ergonomic design providing easy usability
- High capacity Li-ion battery enable long runtime of operation
- Low battery indication
- Powerful motor fills a 25mL pipette in<5 seconds
- 0.45µm replaceable Hydrophobic filter
- Rechargeable during use

Large LCD display showing low battery warning and speed settings

Select the desired coarse speed by press the button +/-

Operating buttons fine-tuned for aspiration and dispensing

Operating button for motor and gravity dispensing

Hydrophobic filter 0.45µm

Silicone adapter for 0.1-100ml pipette tip, fully autoclavable

### Specifications

Aspirate Speeds	25ml < 5s (8 shift)
Dispense Speeds	Motor (8 shift) / Gravity
Battery	Lithium-Ion
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	2~3 Hours
Pipette Types	Glass or Plastic Pipette (0.1-100mL), Pasteur Pipettes
Filter	0.45 µm Hydrophobic
Weight	200g